



BROAD TPLO PLATE 3.5MM			
CAT. NO.	HOLES	SIDE	THICKNESS (MM)
500.941	4	Left	3.7
500.942	4	Right	3.7

Specification

- Broad TPLO plate 3.5mm available in both SS and TT materials.
- The plate contains 4 holes which has 3.7mm thickness and width of 13.5mm. Which is used with Head holes accept 3.5mm locking and 4.0mm cancellous bone screws and shaft holes accept 3.5mm cortex screws.

Indication for use

- These plates are used to treat cranial cruciate ligament (CCL) injuries in dogs.