



DCP BROAD PLATE 3.5MM			
CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.661	7	86	3.6
500.662	8	98	3.6
500.663	9	110	3.6
500.664	10	122	3.6
500.665	11	134	3.6
500.666	12	146	3.6
500.667	13	158	3.6
500.668	14	170	3.6
500.669	16	194	3.6
500.670	18	218	3.6
500.671	20	242	3.6
500.672	22	266	3.6

Specification

- DCP Broad plate available in both SS and TT materials.
- The plate contains 7-22 holes which has 86-266mm in length 3.6mm thickness and width of 12.0mm. Which is used with 3.5mm cortex screws.

Indication for use

- DCP is a metallic plate used in orthopedics for long bone diaphyseal transverse fracture, Iliac body fractures and also used in Intra fractures of ulna and radius
- Fractures of the ischium and pubis

