



H-PLATE 2.0MM			
CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.848	4	6.0	1.0

Specification

- H-Plate 2.0 mm available in both SS and TT materials.
- The plate contains 4 holes which has 6.0 mm in length, 1.0mm thickness and width of 11.0mm.

Which is used with 2.0 mm cortex screws.

Indication for use

- The 2.0 mm Veterinary H-Plate is used for stabilizing small bone fractures in animals, especially in areas requiring low-profile fixation such as the metacarpals, metatarsals, and phalanges. Its H-shaped design provides strong axial and rotational stability for proper bone healing.

