



L-PLATE 2.0MM			
CAT. NO.	HOLES	LENGTH (mm)	Thickness (mm)
500.831	2	17 left	1.0
500.832	2	17 right	1.0
500.833	2	20 left oblique	1.0
500.834	2	20 right oblique	1.0

Specification

- L-Plate 2.0mm available in both SS and TT materials.
 - The plate contains 2 holes which has 17-20mm in length,1.0mm thickness and width of 5.0mm.
- Which is used with 2.0mm cortex screws.

Indication for use