



L-PLATE 2.7 MM			
CAT. NO.	HOLES	LENGTH (mm)	Thickness (mm)
500.838	3	32left	1.0
500.839	3	32right	1.0
500.840	3	34 left oblique	1.0
500.841	3	34 right oblique	1.0

Specification

- L-Plate 2.7mm available in both SS and TT materials.
- The plate contains 3 holes which has 32-34mm in length, 1.0mm thickness and width of 7.0mm. Which is used with 2.7mm cortex screws.

Indication for use