

**LC-DCP PLATE 2.0/1.5MM**

CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.201	4	25	1.2
500.202	5	31	1.2
500.203	6	37	1.2
500.204	7	43	1.2
500.211	6	37	1.5
500.212	7	43	1.5
500.213	8	49	1.5
500.214	9	55	1.5
500.215	10	61	1.5
500.216	12	73	1.5
500.217	14	85	1.5

Specification

- LC-DCP plate available in both SS and TT materials.
- The plate contains 4 to 14 holes which has 25-85 mm in length and 1.2-1.5mm thickness.

Which is used with 1.5mm or 2.0mm cortex screw

Indication for use

- These plates are used for long-bone fractures in small breed dogs and cats.