



LC-DCP PLATE 2.7MM			
CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.241	4	37	2.6
500.242	5	46	2.6
500.243	6	55	2.6
500.244	7	64	2.6
500.245	8	73	2.6
500.246	9	82	2.6
500.247	10	91	2.6
500.248	11	100	2.6
500.249	12	109	2.6
500.250	14	127	2.6
500.251	16	145	2.6

Specification

- LC-DCP plate available in both SS and TT materials.
 - The plate contains 4 to 16 holes which has 37-145mm in length and 2.6mm thickness and width 7.5mm.
- Which is used with 2.7mm cortex screw

Indication for use

- These plates are used for treating fractures, osteotomies, and replantation of small bones, including hand, ankle, and foot, as well as long-bone fractures in the veterinary patients - canines and felines.