



LC-DCP PLATE 3.5 MM

CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.261	2	28	3.3
500.262	3	41	3.3
500.263	4	54	3.3
500.264	5	67	3.3
500.265	6	80	3.3
500.266	7	93	3.3
500.267	8	106	3.3
500.268	9	119	3.3
500.269	10	132	3.3
500.270	11	145	3.3
500.271	12	158	3.3
500.272	14	184	3.3
500.273	15	197	3.3
500.274	16	210	3.3
500.275	18	236	3.3
500.276	20	262	3.3
500.277	22	288	3.3

Specification

- LC-DCP plate available in both SS and TT materials.
- The plate contains 2-22 holes which has 28-288mm in length and 3.5mm thickness and width 11mm.

Which is used with 3.5mm cortex screw

Indication for use

- These plates are used for treating fractures, osteotomies, and replantation of small bones, including hand, ankle, and foot, as well as long-bone fractures in the veterinary patients - canines and felines.