



BROADLCP PLATE 3.5 MM

CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.421	7	94	4.2
500.422	8	107	4.2
500.423	9	120	4.2
500.424	10	133	4.2
500.425	11	146	4.2
500.426	12	159	4.2
500.427	13	172	4.2
500.428	14	185	4.2
500.429	15	198	4.2
500.430	16	211	4.2
500.431	17	224	4.2
500.432	18	237	4.2
500.433	20	263	4.2
500.434	22	289	4.2

Specification

- Board LCP plate available in both SS and TT materials.
 - The plate contains 7-22 holes which has 94-289mm in length and 4.2mm thickness and width 13.5mm.
- Which is used with 3.5mm Locking, 3.5mm Cortex and, 4.0mm Cancellous Bone screws

Indication for use