



LCP BOARD PLATE 5 MM

CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.461	6	107	5.2
500.462	7	125	5.2
500.463	8	143	5.2
500.464	9	161	5.2
500.465	10	179	5.2
500.466	11	197	5.2
500.467	12	215	5.2
500.468	13	233	5.2
500.469	14	251	5.2
500.470	15	269	5.2
500.471	16	287	5.2
500.472	18	323	5.2

Specification

- LCP BOARD PLATE 4available in both SS and TT materials.
- The plate contains 6-18 holes which has 107-323mm in length and 5.2mm thickness and width 17.5mm. Which is used with 4.0mm and 5.0mm Locking, 4.5mm Cortex and, 6.5mm Cancellous Bone screws.

Indication for use