



LCP BOARD PLATE 5.5 MM

CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.481	10	180	6.0
500.482	11	198	6.0
500.483	12	216	6.0
500.484	13	234	6.0
500.485	14	252	6.0
500.486	15	270	6.0
500.487	16	288	6.0
500.488	17	306	6.0
500.489	18	324	6.0

Specification

- LCP BOARD PLATE available in both SS and TT materials.
- The plate contains 10-18 holes which has 180-324mm in length and 6.0mm thickness and width 17.5mm. Which is used with 4.0mm and 5.0mm Locking, 4.5mm and 5.5mm Cortex and, 6.5mm Cancellous Bone screws.

Indication for use