



LCP T - PLATE 3.5 MM 3 HOLES HAED

CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.511	3	50	2.8
500.512	4	57	2.8
500.513	5	67	2.8
500.514	6	77	2.8
500.515	7	87	2.8
500.516	8	97	2.8

Specification

- LCPT-Plate available in both SS and TT materials.
- The plate contains 3-8 holes which has 50-97mm in length and 2.8mm thickness and width 10.00mm. Which is used with 3.5mm Locking, 3.5mm Cortex and, 4.0mm Cancellous Bone screws.

Indication for use

- LCP T-Plate is intended for treatment of physeal fractures of the proximal tibia in foals.