



LOCKING PLATE 2.0MM			
CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.341	4	27	1.2
500.342	5	34	1.2
500.343	6	41	1.2
500.344	7	48	1.2
500.351	6	41	1.5
500.352	7	48	1.5
500.353	8	55	1.5
500.354	9	62	1.5
500.355	10	69	1.5
500.356	12	83	1.5
500.357	14	97	1.5

Specification

- LOCKING PLATE 2.0MM available in both SS and TT materials.
- The plate contains 4-14 holes which has 27-97mm in length and 1.2-1.5mm thickness and width 5.5mm. Which is used with 2.0/1.5mm cortex screw

Indication for use

- These plates are indicated for the treatment of fractures and arthrodesis of small canines and felines.

