



LCP PLATE 2.4 MM			
CAT. NO.	HOLES	LENGTH (MM)	THICKNESS (MM)
500.361	4	32	1.7
500.362	5	40	1.7
500.363	6	48	1.7
500.364	7	56	1.7
500.365	8	64	1.7
500.371	8	64	2.0
500.372	9	72	2.0
500.373	10	80	2.0
500.374	12	96	2.0
500.375	14	112	2.0

Specification

- LCP plate available in both SS and TT materials.
- The plate contains 4-14 holes which has 32-112mm in length and 1.7-2.0mm thickness and width 6.5mm. Which is used with 2.4mm locking, cortex and 2.4mm cortex screws

Indication for use

- These plates are indicated for the treatment of fractures and arthrodesis of small canines and felines.

