



Cat. No.	LCP VETERINARY OSTEOTOMY PLATE	Length MM	Width MM	Thickness MM	Screws	Height
798.920L	LCP VETERINARY OSTEOTOMY PLATE 20°-LEFT	38.5	18.2	3	3.5mm LCP	10
798.925L	LCP VETERINARY OSTEOTOMY PLATE 25°-LEFT	38.5	18.2	3	3.5mm LCP	11.35
798.930L	LCP VETERINARY OSTEOTOMY PLATE 30°-LEFT	38.5	18.2	3	3.5mm LCP	12.75
798.935L	LCP VETERINARY OSTEOTOMY PLATE 35°-LEFT	38.5	18.2	3	3.5mm LCP	13.65
798.920R	LCP VETERINARY OSTEOTOMY PLATE 20°-RIGHT	38.5	18.2	3	3.5mm LCP	10
798.925R	LCP VETERINARY OSTEOTOMY PLATE 25° -RIGHT	38.5	18.2	3	3.5mm LCP	11.35
798.930R	LCP VETERINARY OSTEOTOMY PLATE 30°-RIGHT	38.5	18.2	3	3.5mm LCP	12.75
798.935R	LCP VETERINARY OSTEOTOMY PLATE 35°-RIGHT	38.5	18.2	3	3.5mm LCP	13.65

Specification

- LCP Veterinary osteotomy locking plate Available in both SS & TT materials.
- It is available with 38.5mm length, 18.2mm width, 3mm thickness and 10-13.65mm height. Which is used with 3.5mm Locking screw.
- 20°, 25°, 30° & 35° option

Indication for use

- LCP Veterinary osteotomy locking plate is intended for treatment of coxofemoral joint instability and subluxation in immature dogs prior to onset of osteoarthritis.