

T - PLATE 3.5 MM RIGHT ANGLE



T - PLATE 3.5 MM RIGHT ANGLE			
CAT. NO.	HOLES	LENGTH (mm)	Thickness (mm)
500.810	3	50	1.2
500.811	5	67	1.2
500.812	4	56	1.2
500.813	6	78	1.2
500.814	7	90	1.2
500.815	8	101	1.2
500.816	9	112	1.2

Specification

- **T-Plate 3.5mm right angle** available in both SS and TT materials.
- The plate contains 3-9 holes which has 50-112mm in length, 1.2mm thickness and width of 10.0mm. Which is used with 3.5 cortex screws.

Indication for use