



T-PLATE 2.0 / 1.5 MM, 2 HOLES IN HEAD

CAT. NO.	HAED HOLES	HOLES	LENGTH (mm)	Thickness (mm)
500795	2	4	34	1.5
500.799	2	8	62	1.5

Specification

- T-Plate 2.0/1.5mm available in both SS and TT materials.
- The plate contains 2 head holes and 4-8 holes which has 34-62mm in length,1.5mm thickness and width of 5.0mm. Which is used with 2.0mm and 1.5mm cortex screws.

Indication for use

- These type of plate are used for Radius/Ulna fracture repair in very small dogs.

