


T-PLATE 2.0MM, 9 HOLES



T-PLATE 2.0MM, 9 HOLES

CAT. NO.	HEAD HOLES	HOLES	LENGTH (mm)	Thickness (mm)
500790	3	9	50	1.0
500.791	4	9	50	1.0

Specification

- T-Plate 2.0mm available in both SS and TT materials.
- The plate contains 3-4 head holes and 9 holes which has 50mm in length, 1.0mm thickness and width of 5.0mm. Which is used with 2.0mm cortex screws.

Indication for use

- These type of plate are used for Radius/Ulna fracture repair in very small dogs.