


T-PLATE 3.5 MM / 2.7MM, 2 HOLES IN HEAD



T-PLATE 3.5 MM / 2.7 MM, 2 HOLES IN HEAD				
CAT. NO.	HAED HOLES	HOLES	LENGTH (mm)	Thickness (mm)
500803	2	5	65	2.0
500.805	2	7	85	2.0
500.807	2	9	105	2.0

Specification

- T-Plate 3.5mm/2.7mm available in both SS and TT materials.
- The plate contains 2 head holes and 5-9 holes which has 65-105mm in length, 2.0mm thickness and width of 5.0mm. Which is used with 2.7mm and 3.5 cortex screws.

Indication for use