



TPLO PLATE 2.0MM

CAT. NO.	HOLES	SIDE	Thickness (mm)
500.921	3	Left	2.3
500.922	3	Right	2.3

Specification

- TPLO plate 2.0mm available in both SS and TT materials.
- The plate contains 3 holes which has 2.3mm thickness and width of 5.5mm. Which is used with Head holes 2.0mm locking and 2.0mm cortex screws and shaft holes accept 2.0mm cortex screws.

Indication for use

- Thes type of plate are used for to repair a damaged or torn cranial cruciate ligament (CCL) in the knee of dogs