



TPLO PLATE 2. 7MM

CAT. NO.	HOLES	SIDE	Thickness (mm)
500.931	3	Left	2.5
500.932	3	Right	2.5

Specification

- TPLO plate 2.7mm available in both SS and TT materials.
- The plate contains 3 holes which has 2.5mm thickness and width of 8.0mm. Which is used with Head holes accept 2.7mm locking and 2.7mm cortex screws and shaft holes accept 2.7mm cortex screws.

Indication for use

- Thes type of plate are used for to repair a damaged or torn cranial cruciate ligament (CCL) in the knee of dogs