



TPLO PLATE 3.5MM			
CAT. NO.	HOLES	SIDE	Thickness (mm)
500.936	3	Left	3.7
500.937	3	Right	3.7

Specification

- TPLO plate 3.5mm available in both SS and TT materials.
- The plate contains 3 holes which has 3.7mm thickness and width of 11.0mm. Which is used with Head holes accept 3.5mm locking and 4.0mm cortex screws and shaft holes accept 3.5mm cortex screws.

Indication for use

- Thes type of plate are used for to repair a damaged or torn cranial cruciate ligament (CCL) in the knee of dogs