



CAT. NO. SS	CAT. NO. TT	CORTICAL SCREW 1.5MM (HEXAGONAL HEAD)
100.306	601.206	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 6.0MM
100.308	601.208	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 8.0MM
100.310	601.210	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 10.0MM
100.312	601.212	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 12.0MM
100.314	601.214	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 14.0MM
100.316	601.216	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 16.0MM
100.318	601.218	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 18.0MM
100.320	601.220	CORTICAL SCREW 1.5MM HEXAGONAL HEAD 20.0MM

Specification

- Available in 6 mm to 20 mm length with 1.5mm diameter.
- Available in Stainless steel & Titanium.

Indication for use

- The cortical screw is indicated for use in hard cortical region of the bone. Self-tapping feature eliminates the use of bone tap. These screws are also available with Self-drilling feature.

