


LC - LCP SCREW 2.7 MM (HEXAGONAL HEAD)



CAT. NO. SS	CAT. NO. TT	LC - LCP SCREW 2.7 MM (HEXAGONAL HEAD)
326.010	626.010	LC - LCP SCREW 2.7 MM LENGTH 10 MM (HEXAGONAL HEAD)
326.012	626.012	LC - LCP SCREW 2.7 MM LENGTH 12 MM (HEXAGONAL HEAD)
326.014	626.014	LC - LCP SCREW 2.7 MM LENGTH 14 MM (HEXAGONAL HEAD)
326.016	626.016	LC - LCP SCREW 2.7 MM LENGTH 16 MM (HEXAGONAL HEAD)
326.018	626.018	LC - LCP SCREW 2.7 MM LENGTH 18 MM (HEXAGONAL HEAD)
326.020	626.020	LC - LCP SCREW 2.7 MM LENGTH 20 MM (HEXAGONAL HEAD)
326.022	626.022	LC - LCP SCREW 2.7 MM LENGTH 22 MM (HEXAGONAL HEAD)
326.024	626.024	LC - LCP SCREW 2.7 MM LENGTH 24 MM (HEXAGONAL HEAD)
326.026	626.026	LC - LCP SCREW 2.7 MM LENGTH 26 MM (HEXAGONAL HEAD)
326.028	626.028	LC - LCP SCREW 2.7 MM LENGTH 28 MM (HEXAGONAL HEAD)
326.030	626.030	LC - LCP SCREW 2.7 MM LENGTH 30 MM (HEXAGONAL HEAD)
326.032	626.032	LC - LCP SCREW 2.7 MM LENGTH 32 MM (HEXAGONAL HEAD)
326.034	626.034	LC - LCP SCREW 2.7 MM LENGTH 34 MM (HEXAGONAL HEAD)
326.036	626.036	LC - LCP SCREW 2.7 MM LENGTH 36 MM (HEXAGONAL HEAD)
326.038	626.038	LC - LCP SCREW 2.7 MM LENGTH 38 MM (HEXAGONAL HEAD)
326.040	626.040	LC - LCP SCREW 2.7 MM LENGTH 40 MM (HEXAGONAL HEAD)

Specification

- Conical head facilitates alignment of the locking screw in the threaded plate hole.
- Available in Titanium & Stainless steel.
- Length ranges from 10mm-40mm.

Indication for use

- It is used in the fixation of Small Fragment Wise-Lock Plates.