

LC - LCP SCREW 5.0 MM (HEXAGONAL HEAD)



CAT. NO. SS	CAT. NO. TT	LC - LCP SCREW 5.0 MM LENGTH HEXAGONALHEAD
328.510	628.510	LC-LCP Screws Ø 5.0 mm x 10mm
328.512	628.512	LC-LCP Screws Ø 5.0 mm x 12mm
328.514	628.514	LC-LCP Screws Ø 5.0 mm x 14mm
328.516	628.516	LC-LCP Screws Ø 5.0 mm x 16mm
328.518	628.518	LC-LCP Screws Ø 5.0 mm x 18mm
328.520	628.520	LC-LCP Screws Ø 5.0 mm x 20mm
328.522	628.522	LC-LCP Screws Ø 5.0 mm x 22mm
328.524	628.524	LC-LCP Screws Ø 5.0 mm x 24mm
328.526	628.526	LC-LCP Screws Ø 5.0 mm x 26mm
328.528	628.528	LC-LCP Screws Ø 5.0 mm x 28mm
328.530	628.530	LC-LCP Screws Ø 5.0 mm x 30mm
328.532	628.532	LC-LCP Screws Ø 5.0 mm x 32mm
328.534	628.534	LC-LCP Screws Ø 5.0 mm x 34mm
328.535	628.535	LC-LCP Screws Ø 5.0 mm x 35mm
328.536	628.536	LC-LCP Screws Ø 5.0mm x 36mm
328.538	628.538	LC-LCP Screws Ø 5.0 mm x 38mm
328.540	628.540	LC-LCP Screws Ø 5.0 mm x 40mm
328.542	628.542	LC-LCP Screws Ø 5.0 mm x 42mm
328.544	628.544	LC-LCP Screws Ø 5.0 mm x 44mm
328.545	628.545	LC-LCP Screws Ø 5.0 mm x 45mm
328.546	628.546	LC-LCP Screws Ø 5.0 mm x 46mm
328.548	628.548	LC-LCP Screws Ø 5.0 mm x 48mm
328.550	628.550	LC-LCP Screws Ø 5.0 mm x 50mm
328.552	628.552	LC-LCP Screws Ø 5.0 mm x 52mm
328.554	628.554	LC-LCP Screws Ø 5.0 mm x 54mm



328.555	628.555	LC-LCP Screws Ø 5.0 mm x 55mm
328.556	628.556	LC-LCP Screws Ø 5.0 mm x 56mm
328.558	628.558	LC-LCP Screws Ø 5.0 mm x 58mm
328.560	628.560	LC-LCP Screws Ø 5.0 mm x 60mm
328.565	628.5645	LC-LCP Screws Ø 5.0 mm x 65mm
328.570	628.570	LC-LCP Screws Ø 5.0 mm x 70mm
328.575	628.575	LC-LCP Screws Ø 5.0 mm x 75mm
328.576	628.576	LC-LCP Screws Ø 5.0 mm x 76mm
328.580	628.580	LC-LCP Screws Ø 5.0 mm x 80mm
328.582	628.582	LC-LCP Screws Ø 5.0 mm x 82mm
328.584	628.584	LC-LCP Screws Ø 5.0 mm x 84mm
328.585	628.585	LC-LCP Screws Ø 5.0 mm x 85mm
328.590	628.590	LC-LCP Screws Ø 5.0 mm x 90mm
328.600	628.600	LC-LCP Screws Ø 5.0 mm x 100mm

Specification

- Available in 10mm to 100mm length with 5.0mm diameter.
- Fully threaded shaft.
- Available in both Stainless Steel & Titanium.

Indication for use

- It is used in the fixation of Small Fragment Wise-Lock Plates.

