

LATERAL MASS / POSTERIOR CERVICAL POLY SCREWS 3.5 MM



SR. NO.	CAT. NO. TT	METAL	LATERAL MASS / POSTERIOR CERVICAL POLY SCREWS 3.5 mm
1	692.510	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 10 mm
2	692.512	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 12 mm
3	692.514	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 14 mm
4	692.516	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 16 mm
5	692.518	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 18 mm
6	692.520	TT	Lateral Mass / Posterior cervical Polyscrews 3.5 mm / 20 mm
7	692.522	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 22 mm
8	692.524	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 24 mm
9	692.526	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 26 mm
10	692.528	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 28 mm
11	692.530	TT	Lateral Mass / Posterior cervical Poly screws 3.5 mm / 30 mm

Specification

- Lateral Mass / Posterior cervical Poly screws Available in Titanium.
- Available length from 10,12,14,16,18,20,22,24,26,28 & 30 with 3.5mm diameter.

Indication for use

- Lateral Mass Poly axial screw is used for stabilization of the occipitocervical junction.

