



CAT. NO.	METAL	Gouge
1761.525	SS	Gouge with Fibre Handle - 25 mm, Straight
413.050	SS	Gouge with Fibre Handle - 5 mm, Straight
413.051	SS	Gouge with Fibre Handle - 5 mm, Curved
413.060	SS	Gouge with Fibre Handle - 6 mm, Straight
413.061	SS	Gouge with Fibre Handle - 6 mm, Curved
413.100	SS	Gouge with Fibre Handle - 10 mm, Straight
413.101	SS	Gouge with Fibre Handle - 10 mm, Curved
413.120	SS	Gouge with Fibre Handle - 12 mm, Straight
413.150	SS	Gouge with Fibre Handle - 15 mm, Straight
413.151	SS	Gouge with Fibre Handle - 15 mm, Curved
413.160	SS	Gouge with Fibre Handle - 16 mm, Straight
413.161	SS	Gouge with Fibre Handle - 16 mm, Curved
413.180	SS	Gouge with Fibre Handle - 18 mm, Straight
413.181	SS	Gouge with Fibre Handle - 18 mm, Curved
413.201	SS	Gouge with Fibre Handle - 20 mm, Curved
413.202	SS	Gouge with Fibre Handle - 20 mm, Straight
413.24 0	SS	Gouge with Fibre Handle - 24 mm, Straight
413.241	SS	Gouge with Fibre Handle - 24 mm, Curved
413.251	SS	Gouge with Fibre Handle - 25 mm, Curved
413.300	SS	Gouge with Fibre Handle - 30 mm, Straight
413.301	SS	Gouge with Fibre Handle - 30 mm, Curved

Specification

- Gouge is available in Stainless steel.
- It is available in 5,6,10,12,15,16,18,20,24,25 & 30mm length.

Indication for use

- A bone gouge is a hollow chisel that is used to cut into and remove portions of bone. They can be used in notchplasty, to demarcate and remove notch fragments, and in other orthopedic procedures to contour bone or remove cement.

